

# SPECIFICATIONS

## POROLUX™ Pressure Step Stability Series



### PRODUCT OVERVIEW

### POROLIQT™ 1000AQ

### POROLIQT™ 1000ML

Measuring principle	Pressure driven	Pressure driven
Max. Pressure	40 bar	40 bar
Measurement range	0.3 µm - 2 nm	1 µm - 2 nm
Flow range	0.16 µl - 10 ml/min	1 µl - 10 ml/min
Sample holders	Flat sheet 25mm diameter Hollow fiber	Flat sheet 25mm diameter Hollow fiber
Special sample holders	Available upon request	Available upon request
Pressure sensors	50 bar	50 bar
Displacement liquid	Water	Multiple liquids