

SPECIFICATIONS

SentinelPro™



PHYSICAL

Height	38.1 cm (15 in.)
Width	25.4 cm (10 in.)
Depth	63.5 cm (25 in.)
Weight	13.2 kg (29lbs.)
System volume	PI SENTINEL Shape Module allows for different sample volumes.

ELECTRICAL

Voltage	100 - 240 VAC.
Frequency	50 - 60 Hz.

ENVIRONMENT

Temperature	10 to 45 °C for operation. -10 to 55 °C for storing or shipping.
Humidity	20 to 80% relative, without condensation.

COMPUTER REQUIREMENTS

If not using the supplied All-in-One computer

Processor	2 GHz CPU speed or superior.
Memory	1 GB of RAM (minimum).
Hard Dish Space	60 GB hard drive.
Monitor	24" LED Touchscreen
Media Drive	Writable CD or DVD drive.
USB interface	2 USB, 3 connectors.

Due to continuous improvements, specifications are subject to change without notice.